

# TRUGENE HBV Genotyping Assay

Identify Viral Genotype and Interrogate Hepatitis B Viral Mutations

Answers for life.

**SIEMENS**

For Research Use Only. Not For Use  
In Diagnostic Procedures.





## Introduction

The Hepatitis B viral genome is a partially double-stranded DNA molecule, approximately 3200 bases in length, packed inside a nucleocapsid structure composed of the viral core (C) protein (Figure 1). This nucleocapsid is surrounded by the viral envelope, comprising of a mixture of lipids and viral structural proteins or surface antigens (HBsAg, Pre-S1 and Pre-S2).

Virus replication is carried out through enzymatic activity of the viral polymerase or P protein.

This multifunctional protein has domains possessing DNA synthesis priming, reverse transcriptase (RT)/DNA polymerase and ribonuclease (RnaseH) activities<sup>1</sup>. Within the RT region, the protein can be further subdivided into several functional domains (A – F), based on analogy with other viral polymerases. The organization of the HBV genome is such that the same region that codes for the viral P protein also codes for the viral structural proteins, but in an alternate translational reading frame<sup>2</sup> (Figure 2).

# TRUGENE® HBV Module 2.0.

## Identify and Interrogate Hepatitis B Viral Mutations

### Antiviral Drug Development

The P protein and the reverse transcriptase/DNA polymerase domain are the major functional units for viral replication, and as such have been targets in HBV antiviral drug development. A number of new therapeutic agents that inhibit viral replication are currently approved and in clinical trials.

Given the resistance issues identified with lamivudine therapy and the beneficial experience of using combinations of several drugs in the HIV setting, it is likely that new agents will be used in combination with lamivudine to attempt to suppress the development of resistance.

Resistance to lamivudine is a consequence of the low fidelity of the viral polymerase, the high viral load, high viral replicative capacity and the long course of therapy.<sup>10,11,12</sup> The incidence of drug resistance increases with length of therapy, so that typically a significant number of patients will have resistant virus after a year, rising after three years.

### Molecular Variants of HBV

Originally, there were two different numbering schemes for assigning mutations within the reverse transcriptase (RT) domain of HBV. These schemes are known in the scientific literature as the genotype-dependent numbering and the genotype-independent numbering scheme.

The lack of a universally accepted numbering convention led to inconsistencies in the published scientific literature. As a result, the genotype-independent scheme has emerged as a method to standardize the way in which HBV mutations in the RT domain are numbered.<sup>3</sup>

The TRUGENE HBV Module 2.0 has adopted that widely accepted nomenclature. This classification system is based on analysis of the complete genomic sequence of the virus. Eight discrete genotypes (A-H) have been identified.<sup>4,5,6,7</sup>

It has been suggested that the viral genotype may correlate with differences in clinical features of HBV infection.<sup>8,9</sup>

The region of the genome that is used to assign viral genotype is also the region that is associated with the development of drug resistance.




Figure 1. The Hepatitis B Virus.




Figure 2. The HBV Genome.

Reverse Transcriptase Mutations Identified by the TRUGENE HBV Genotyping Assay
A181T
A181IV + M204I
A181V/T
I169T + V173L + L180M + <b>T184G</b> + <b>S202I</b> + M204V
I169T + V173L + L180M + M204V + <b>M250V</b>
I169 T + V173L + L180M + M204V
L180M + M204V/I/S
L180M + A194T + M204V
L180M + M204V/I
L180V/I + M204I [non-A genotype]
M204I
N236T
Q215S + L180M + M204V
T184S + L180M+ M204V
V173L + L180M + M204V
V214A / Q215S
Various combinations of mutations at codons <b>184, 202</b> and <b>250</b>

List is not comprehensive. Any mutation that occurs within the primer region (Rt99-rt280 and (s101-s237) will be identified and listed in one of the four Mutation Profile columns.



The OpenGene DNA Sequencing System.

## TRUGENE HBV Genotyping Assay – Two Tests in One Kit

The TRUGENE HBV Assay from Siemens Healthcare Diagnostics allows researchers to identify the HBV genotype and HBV mutations at the same time using direct sequence analysis. This improves workflow efficiency and reduces overall costs normally associated with separate HBV genotyping and mutational analysis. The assay runs on Siemens' OpenGene® DNA Sequencing System, offering a sequence-based solution for determination of:

- The genotype
- Viral mutations, highlighting those that differ from a consensus sequence and that current scientific literature has identified

This Research Use Only assay amplifies the viral genome directly from plasma or serum samples and sequences the region coding for the viral RT gene and the central portion of HBsAg. Amplified material is then added to a set of four CLIP™ sequencing reaction tubes and labeled via the CLIP sequencing reaction (Figure 3).




Figure 3. Timeline for TRUGENE HBV genotyping.

# A fully integrated HBV Genotyping Solution: Hardware, Chemistry, Software and Report

Traditional approaches to HBV genotyping have involved "home brew" sequence-based phylogenetic analysis<sup>11</sup>, differential hybridization<sup>12</sup>, fragment length polymorphism<sup>13</sup>, or subjective hybridization based assays. The TRUGENE HBV Genotyping Assay is a standardized and validated kit allowing researchers to directly sequence and interrogate the virus for a more comprehensive viral analysis. This provides the laboratory more confidence in the data generated for HBV genotyping and mutational analysis.

## CLIP Sequencing

CLIP is a proprietary sequence chemistry that produces bi-directional sequences using two fluorescently-labeled DNA primers which increases the overall confidence of the base calls. The sequencing reaction is initiated with the addition of the sample and a thermostable DNA polymerase. As the reaction mixture proceeds through each thermal cycle, primers hybridize to template DNA and are extended, then terminated along the target DNA sequence allowing bi-directional sequence in the same tube.

## Analysis

The OpenGene DNA Sequencing System acquires sequence data in real time, and each pair of forward and reverse sequences are combined and aligned with stored reference sequences. The TRUGENE HBV Genotyping Kit Includes reagents to PCR amplify a portion of the HBsAg (s101-s237) and overlapping polymerase protein (rt99-rt280). Software assigns the viral genotype and the mutations and polymorphisms present in the sample.

The TRUGENE HBV software module contains sequences that correspond to the surface antigen and polymerase regions for genotype A through H and a universal mutation reporting reference sequence for comparison.

## Reporting

The TRUGENE HBV Genotyping Report shows the closest consensus sequence allowing the determination of the viral genotype and a mutation table showing the published mutations from key scientific journal articles.

## Conclusion

The TRUGENE HBV Assay from Siemens Healthcare Diagnostics runs on the OpenGene DNA Sequencing System and uses sequencing technology that overcomes the traditional issues faced when performing DNA sequencing.

The assay provides a determination of HBV viral genotype and interrogates known mutations – essentially providing two answers from one test.



**TRUGENE® HBV Genotyping Report**

Sample ID: EN-7209      Laboratory: Your laboratory Name  
 Date Drawn: Aug 17, 2005      4321 Elm Street  
 Case ID: 0001      Anytown, USA 50535  
 Case Name: Herb Nichols  
 Investigator: Dr. Jane Doe  
 Institution: City Hospital  
 Reported by: Anna Lyst  
 Report Date: 2009/04/14      Tel: 800-434-2447 Fax: 510-705-5902

---

**HBV Genotype: A**  
 (i.e. the top match w. \_A\_302763)

Type	Name	% Homology	Score / 133
A	X02763	99.0	132
A	X51970	97.8	128
A	AF090842	95.7	115
D	X65259	94.5	109
B	D00329	90.1	102
B	AF100309	92.7	100
G	AF180501	91.7	95
G	AB043410	91.7	95
G	AF405706	91.6	94
B	AB033554	91.2	92

The homology is the percentage of identical bases within the sequenced region, between the sample and the isolate. The score indicates the number of identical bases at known polymorphic positions. The value in the % Homology column is determined from the entire region sequenced. The value in the Score column is determined from the Surface Antigen gene only.

Authorized Reviewer: \_\_\_\_\_ Signature: \_\_\_\_\_ Date: \_\_\_\_\_  
 Name: \_\_\_\_\_ Title: \_\_\_\_\_

HBV Genotype analysis by DNA sequencing  
 For Research Use Only. Not for use in diagnostic procedures.  
 © 2008 Siemens Healthcare Diagnostics Inc. All Rights Reserved.  
 GeneOpGen(TM) 4.1 (2008/05/28) / HBV\_S\_RT\_3002.gr1 Fri Feb 27 2009 10:00 AM Page 1 of 3

**TRUGENE® HBV Genotyping Report**

Sample ID: EN-7209      Laboratory: Your laboratory Name  
 Date Drawn: Aug 17, 2005      4321 Elm Street  
 Case ID: 0001      Anytown, USA 50535  
 Case Name: Herb Nichols  
 Investigator: Dr. Jane Doe  
 Institution: City Hospital  
 Reported by: Anna Lyst  
 Report Date: 2009/04/14      Tel: 800-434-2447 Fax: 510-705-5902

---

**Mutation Profile\***

Published Mutations**	Silent Mutations (at all positions)	Coding changes not at Published Sites**	Unpublished Coding Changes at Published Sites**	
None	87 386 (GAG) D1902 (CAA) 400 (GAT) D1340 (GAC) 465 (CTG) L187 (CTG) 468 (TTG) L187 (CTG) 502 (AAC) N189 (AAC) 605 (ATT) H188 (ATA) 688 (CTG) L188 (CTG) 683 (CAA) D2192 (CAA) 881 (TTG) L216 (CTG) 783 (TTG) L238 (TTA) 713 (ATT) N209 (AAC) 718 (ACC) T247 (ACA) 720 (TAC) R249 (AAG) 730 (GCG) G240 (GGT) 733 (TAC) T249 (TAT) 738 (CTT) L247 (CTA) 739 (GAA) G249 (GAT) 754 (TAT) V257 (TAC) 772 (GTT) S269 (GAA) 778 (TTA) L268 (TTG) 781 (CTC) P257 (CCA) 802 (AAN) K266 (AAG) 805 (AAC) G270 (AAG) 817 (TTT) P270 (TTC) 822 (GAA) R270 (AAG) 828 (TTG) L270 (CTT) 832 (CTA) N270 (GTT) 835 (AA) T297 (AA)	331 (TAT) Y1119F (TAT) 334 (CTT) Y1120F (TCT) 362 (ACT) T1189 (ACT) 384 (AAC) N2101 (ACT) 384 (TAC) Y1200 (AAC) 378 (TAC) Y1200 (AAC) 388 (GAG) N199 (GAC) 413 (AAC) N1993 (CAA) 481 (TTG) L199 (TTG) 481 (TTG) F197 (TAC) 481 (CTC) H198 (CTC) 481 (CTT) L199 (GCT) 481 (TTT) P211 (TAT) 681 (CTT) S234 (GCT) 676 (ATT) S244 (GTT) 781 (GAA) E262 (GAT) 805 (CTC) L269 (ATC) 811 (CAA) G271 (CAC)	5 367 (ATG) M1088L (ATG) 328 (CTT) L1102 (ATT) 337 (ACT) T1130 (ACT) 340 (TCA) S1147 (ACA) 340 (AAC) N1070 (AAC) 376 (ATT) H187 (ACT) 376 (ACT) Y1170 (ACT) 427 (TCC) S1437 (ACC) 479 (AAC) H186 (AAG) 481 (TTG) P117 (TAC) 586 (AAC) N1870 (AAC) 605 (TTG) L209 (GTT) 607 (TTA) L210 (ATA)	None

HBV Genotype analysis by DNA sequencing  
 For Research Use Only. Not for use in diagnostic procedures.  
 © 2008 Siemens Healthcare Diagnostics Inc. All Rights Reserved.  
 GeneOpGen(TM) 4.1 (2008/05/28) / HBV\_S\_RT\_3002.gr1 Fri Feb 27 2009 10:00 AM Page 2 of 3

TRUGENE HBV reporting.

# OpenGene DNA Sequencing System

Our OpenGene DNA Sequencing System runs a panel of TRUGENE Genotyping Assays to facilitate semi-automated DNA sequencing and/or genotyping.

- Backward compatibility of data for easy conversion from OpenStep (UNIX-based system)

- Fully integrated, dynamic genotype determination
- Accessible – allows research laboratories to initiate DNA sequencing operations with minimal equipment and a small footprint
- Scalable – accommodates up to eight Long-Read Towers per MAC OS<sup>®</sup>X unit for easy scale-up
- Quality – CLIP sequencing chemistry is optimized to improve workflow by eliminating several steps associated with other sequencing methodologies
- MAC OSX operating system and OpenGene software provide enhanced data management and easy networking



TRUGENE HBV Genotyping Assay For Research Use Only.  
 Not For Use In Diagnostic Procedures.

## References:

1. Maramorosch, K., et al. (1999). Advances in Virus Research; published by Academic Press, pp 95-96.
2. Tillmann, H. L., et al. (1999). "Mutational pattern of hepatitis B virus on sequential therapy with famciclovir and lamivudine in patients with hepatitis B virus reinfection occurring under HBIg immunoglobulin after liver transplantation." Hepatology; 30(1):244-256.
3. Stuyver, L. J., et al. (2001). "International Committee Nomenclature for antiviral-resistant human hepatitis B virus mutations in the polymerase region." Hepatology 33(3):751-757.
4. Norder, H., B. Hammas, et al (1990). "Typing of hepatitis B virus genomes by a simplified polymerase chain reaction." J Med Virol 31(3): 215-21.
5. Stuyver, L., S. De Gendt, et al. (2000). "A new genotype of hepatitis B virus: complete genome and phylogenetic relatedness." J Gen Virol 81 Pt 1: 67-74.
6. Norder, H., Courouze A. M., et al. (1994). "Complete genomes, phylogenetic relatedness and structural proteins of six strains of the hepatitis B virus, four of which represent new genotypes." Virology 198: 489-503.
7. Naito, H., et al. (2001). "Rapid and specific genotyping system for hepatitis B virus corresponding to six major genotypes by PCR using type-specific primers." J Clin Microbiol 39: 362-4.
8. Kao, J, H., P. J. Chen, et al. (2000). "Hepatitis B genotypes correlate with clinical outcomes in patients with chronic hepatitis B." Gastroenterology 118(3): 554-9.
9. Fujie, H., K. Moriya, et al. (2001). "Hepatitis B virus genotypes and hepatocellular carcinoma in Japan." Gastroenterology 120(6): 1564-5.
10. Dienstag, J. L., Schiff, E. R., et al (1999). "Lamivudine as initial treatment for chronic hepatitis B in the United States." New England Journal of Medicine 341: 1256-1263.
11. Liaw, Y. F., R. N. Chien, et al. (1999). "Acute exacerbation and hepatitis B virus clearance after emergence of YMDD motif mutation during Lamivudine therapy." Hepatology 30(2): 567-72.
12. Leung, N. W. Y., et al. (2001). "Extended Lamivudine treatment in patients with chronic hepatitis B enhances hepatitis B e antigen seroconversion rates: Results after 3 years of therapy." Hepatology 33(6): 1527-32.

Siemens Healthcare Diagnostics, the leading clinical diagnostics company, is committed to providing clinicians with the vital information they need for the accurate diagnosis, treatment and monitoring of patients. Our comprehensive portfolio of performance-driven systems, unmatched menu offering and IT solutions, in conjunction with highly responsive service, is designed to streamline workflow, enhance operational efficiency and support improved patient care.

CLIP, MicroCel, OpenGene, TRUGENE and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc. All other trademarks and brands are the property of their respective owners.

Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

### Siemens Global Headquarters

Siemens AG  
Wittelsbacherplatz 2  
80333 Muenchen  
Germany

### Global Siemens Healthcare Headquarters

Siemens AG  
Healthcare Sector  
Henkestrasse 127  
91052 Erlangen  
Germany  
Telephone: +49 9131 84 - 0  
[www.siemens.com/healthcare](http://www.siemens.com/healthcare)

### Global Division

Siemens Healthcare Diagnostics Inc.  
1717 Deerfield Road  
Deerfield, IL 60015-0778  
USA  
[www.siemens.com/diagnostics](http://www.siemens.com/diagnostics)

[www.siemens.com/diagnostics](http://www.siemens.com/diagnostics)